

### Comfort - Functionality - Durability

#### Features:

- The Lumex FR566G Deluxe Clinical Care Recliner features three positions for patient comfort. The fourth position, a caregiver-activated Trendelenburg (or Shock) Position, is easily actuated from either side of the recliner by a conveniently located foot pedal
- The wide, contoured seat and back, along with an integrated footrest, provide an enhanced level of patient comfort
- Wide range of fabrics / colors
- Quick Ship available
- Wall saver armrest caps increase durability of the armrests and protect walls from marring
- Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with one directional locking caster

- Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Side table surface is smooth and features tapered edges for easy cleaning
- 350 lb maximum weight capacity, evenly distributed

#### Options:

- Heat and Massage
- Fully Upholstered Armrests
- IV Pole Holder / IV Pole
- Headrest Cover
- Universal Pillow
- Available with Drop Arm Feature (Model FR566DG)

#### Dimensions

Overall Height ..... 50" / 127 cm  
 Overall Width ..... 31" / 79 cm  
 Overall Depth ..... 43" / 109 cm  
 Width Between Arms ..... 21" / 53 cm  
 Arm Height From Seat ..... 8" / 20 cm  
 Seat Height ..... 22" / 56 cm  
 Seat Width ..... 20.5" / 52 cm

Seat Depth ..... 20" / 51 cm  
 Depth Fully Reclined  
 In Trendelenburg ..... 75" / 191 cm  
 Width With Both Side Tables  
 in Upright Position ..... 51" / 130 cm  
 Maximum Weight Capacity  
 Evenly Distributed ..... 350 lb / 158.8 kg



*Upholstered inside panels*

*Smooth tables with no-seam tapered edge*



*Easily removed side panels for cleaning and maintenance*



*Trendelenburg Position can be actuated from either side of the recliner by a convenient foot pedal*



*Multiple Positions for Enhanced Patient Comfort*