

## EXAM ROOM TABLES



All 4200 Series exam tables have metal pass-through drawers rather than the typical plastic, so items can be accessed on both sides.





continuous lock, the most in the industry. A procedure tray is independent of the leg rest

and may be accessed easily at any time. An all-metal construction for the base table provides longevity and recyclability.

- Three Pass Through Drawers
- Two Storage Drawers
- Pull-Out Leg Rest with Integrated Treatment Tray
- Five Position Continuous Lock Stirrups
- Powder Coat Finish Foghorn
- Pelvic Tilt

Exam Table with Three Pass-Through Drawers and Two Storage Drawers and Stirrups

420	Λ \	~~	C ~ ~	: ~ ~
471			Ser	IIES
	,		001	

Height: 32 in / 81 cm

Overall Length: 72 in / 183 cm (Leg Rest Extended)

Width: 27 in / 69 cm

Foot Step: W 18 in / 46 cm X L 17 in / 43 cm

Paper Roll Holder Widths:

18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Drawer Size: W 18 in / 46 cm X L 16.5 in / 42 cm X D 4 in / 10 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 330 lbs / 150 kg

Replacement Top: GF4200RT-XX

Options:

4200-E0-XX **Electrical Outlet** 

4200-DW-XX Drawer Warmer (Includes Electrical Outlet)

4200-TB133-XX TB133 Vinyl

- Dual Pedestal Base with DC Electromechanical Drive
- Pull-Out Leg Rest with Integrated Treatment Tray
- Power Assisted Back
- Five Position Continuous Lock Stirrups
- Hand Pendant: Hi-Low, Back, Flat, Return to Chair
- Powder Coat Finish Foghorn
- Pelvic Tilt
- **UL Listed**

High/Low Power Table with Power Assisted Back, Hand Control and Stirrups

## 4600-XX Series

Height: 18.5 in / 47 cm to 35.5 in / 90 cm\*

Overall Length: 72 in / 183 cm (Leg Rest Extended)

Width: 27 in / 69 cm

Leg Rest: W 18 in / 46 cm X L 17 in / 43 cm

Paper Roll Holder Widths:

18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 700 lbs / 317 kg (EVENLY DISTRIBUTED)

Shipping Weight: 395 lbs / 179 kg

## Replacement Top: GF4600RT-XX

Options:

4600-FS-XX Foot Control 4600-PS-XX Physician Storage 4600-AA-XX Assist Arms

4600-TB-XX TB133 Vinyl 4600-DW-XX Drawer Warmer

4600-FL-XX All above options except TB133

## Standard Color Options(-xx)

Available in 24 color choices to complement most decors

WH	DO	PM	SH	AR	TE
Super White	Doe	Palm	Shadow	Artichoke	Teal
NO	DB	PU	WW	BT	GR
Northwoods	Dusty Blue	Plum	Wedgewood	British Tan	Gray
AL	ST	PL	TL	MO	OG
Alabaster	Saddle Tan	Pearl	Tea Leaf	Moss	Ocean Gray
AQ	BU	VI	NA	WG	BL
Agua Green	Blue	Violet	Navy	Warm Gray	Black

Note: Due to the many variations in monitors and browsers, color samples may appear different. Swatch Card is available from your representative.











\*+- .5 in