

TREATMENT TABLES









Special Features

Our pass-through drawers can be opened from either side of the table, so you'll always have easy access to supplies no matter where your table sits.



Treatment Table, Flat Top

4100-XX Series

Finish: Powder Coat, Foghorn

Overall Length: 76 in / 193 cm

Width: 28 in / 71 cm

Height: 32 in / 81 cm

Paper Roll Holder Widths:

18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 170 lbs / 77 kg

Treatment Table, Flat Top, Shelf and Two Drawers

4103-DD-XX Series

Finish: Powder Coat, Foghorn

Overall Length: 76 in / 193 cm

Width: 28 in / 71 cm

Height: 32 in / 81 cm

Paper Roll Holder Widths:

18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 170 lbs / 77 kg

Treatment Table with Adjustable Back and Contour Top

4105-XX Series

Finish: Powder Coat, Foghorn

Overall Length: 76 in / 193 cm

Width: 28 in / 71 cm

Height: 32 in / 81 cm

Paper Roll Holder Widths:

18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 190 lbs / 86 kg

Treatment Table with Adjustable Back, Contour Top, Shelf and Two Drawers

4105-DD-XX Series

Finish: Powder Coat, Foghorn

Overall Length: 76 in / 193 cm

Width: 28 in / 71 cm

Height: 32 in / 81 cm

Paper Roll Holder Widths:

18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 210 lbs / 95 kg

Standard Color Options(-xx)

Available in 24 color choices to complement most decors

Available III 24 color choices to complement most decors.					
WH	DO	PM	SH	AR	TE
			-		
Super White	Doe	Palm	Shadow	Artichoke	Teal
NO	DB	PU	WW	BT	GR
NU	νD	PU	VVVV	DI	un
Northwoods	Dusty Blue	Plum	Wedgewood	British Tan	Gray
Λ.Ι.	ОТ	DI	TI	1.40	00
	ST	PL	IL	MO	OG
	Saddle Tan	Pearl	Tea Leaf	Moss	Ocean Gray
4.0	511	. "		1410	5.
AQ	BU	\bigvee	NA	WG	BL
Agua Green	Blue	Violet	Navy	Warm Gray	Black

Note: Due to the many variations in monitors and browsers, color samples may appear different. Swatch Card is available from your representative.







