



## Platinum Collection Bath Seat with Arms - Retail

- · Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- · Available with or without backrest
- · Available with or without support arms
- · Tool Free design
- · Designed to meet the requirements of HCPCS code: E0245
- Maximum Weight Capacity 350 lb EVENLY DISTRIBUTED
- · Limited lifetime warranty

7921RA with backrest 3/cs 7931RA without backrest 4/cs

## Specifications

Overall Width	
Overall Depth	18"
Seat Width	20"
Seat Depth	
Minimum Seat Height	
Maximum Seat Height	20"
Width Between Arms	
Maximum Weight Capacity	350 lb EVENLY DISTRIBUTED







