

Two Button Folding Walker with 5" Wheels

Features:

- Easy push dual-release folding mechanism enables user to fold walker sides independently
- Sturdy 1" diameter steel construction ensures maximum strength
- PVC contoured hand grip
- Comes with 5" wheels (Replacement leg set with wheels item #716280-LEGS)
- Rear glide tips allow for walker to move easily and smoothly over most surfaces
- Maximum Weight Capacity
300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0143
- One (1) year limited warranty



Item	Description	UOM
716280A-4	Two Button Folding Walker	4 ea/cs
716280-LEGS	Replacement Leg Set	1 set
716280-TIPS	Replacement Tips	1 pr

Specifications

Width Inside Handgrip.....	17.5" (44.45cm)
Overall Width.....	24.5" (62.23cm)
Depth Front to Rear.....	20.25" (51.44cm)
Minimum Handle Height.....	32.5" (82.55cm)
Maximum Handle Height.....	39.5" (100.33cm)
Approx User Height.....	5'4" – 6'2"
Weight Capacity	300 lb (136 kg)
	Evenly Distributed
Product Weight.....	7.5 lb (3.45 kg)