

T.E.N.S. (Transcutaneous Electrical Nerve Stimulation) is indicated under a physician's prescription for the symptomatic relief and management of chronic (long-term) pain and the treatment of post-operative or post-traumatic pain.



Standard Features:

- Dual channels
- Three modes:
Standard, Burst, Modulation
- Includes: TENS unit, soft carrying case, instructions and accessories (four electrodes and lead wires)
- Also required (not included):
One 9-Volt alkaline battery
- Additional Timer Feature: GF-3T

Federal Law (USA) restricts this device to sale by or on the order of a physician.



Item	Description	UOM
GF-3	Analog TENS	1 ea
GF-3T	Analog TENS with Timer	1 ea
9260	2"x 2" Replacement Electrodes	40/bx
9261	Replacement Butterfly Electrodes	1 ea
GF-45	Replacement Lead Wires	1 ea

