

Everest & Jennings introduces the Traveler[®] L3 Plus as the lightweight wheelchair addition to the Traveler[®] line, a brand respected for performance and durability. The Traveler[®] L3 Plus lightweight frame weighs less than 36 pounds making it the ideal K0003 Medicare coded wheelchair built for short and long term use as well as rental needs.

This lightweight wheelchair comes standard with non-marking polyurethane tires, industry compatible front rigging and adjustable anti-tippers along with a detachable 2" seat extension with hardware. It features a dual position axle, multi-position caster forks for 18 or 20 inch, 16.3" or 17 1/4" seat height adjustment and offers flip-back desk, adjustable height and full arm options for easy lateral movement. The Traveler[®] L3 Plus is manufactured with the same commitment to quality Everest & Jennings is known for.



2" seat extension with hardware comes standard



Anti-tippers come standard

Standard Features

- Weighs less than 36 lb without front rigging
- Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary
- Two position axle and caster fork provides 18" to 20" seat height in regular models and seat height 16.3" to 17.9" in the XVI models
- Hammertone finish and carbon steel frame provide strength and durability
- Black nylon upholstery with standard chart pocket on back
- Includes fastener tab under seat connecting both seat and back upholstery to eliminate gap
- Anti-tippers come standard
- Non-marking polyurethane tires with 24" rear wheels with 8" casters in the regular models and 20" rear wheels with 6" casters in the XVI models
- Tool-free adjustable front rigging with heel loops standard on swingaway footrests
- Flip-back padded arms available in fixed desk length, adjustable height desk and detachable flip-back full length
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003



Warranty

Five (5) year limited warranty



Chart pocket / Traveler[®] L3 Plus XVI model













Full length arms



Adjustable height arms

**This product is not made with natural rubber latex*

Traveler® L3 Plus		Dimensions and Weights									
		Seat			Back	Overall				Weight	
		Width	Depth	Height	Height	Depth	Height	Width	Folded Width	Without front rigging	Max Capacity
											
3F012120	Flip-Back Desk Arm, Swingaway detachable footrest	18"	16"*	18" or 20"	17 ¼"	29"	35½"	26"	12 ¼"	34 lb	300 lb
3F012130	Flip-Back Desk Arm, Elevating legrest	18"	16"*	18" or 20"	17 ¼"	29"	35½"	26"	12 ¼"	34 lb	300 lb
3F012140	Flip-back Full Arm, Swingaway Footrest	18"	16"*	18" or 20"	17 ¼"	29"	35 ½"	26"	12 ¼"	34 lb	300 lb
3F012220	Flip-Back Desk Arm, Swingaway detachable footrest	16"	16"*	18" or 20"	17 ¼"	29"	35½"	24"	12 ¼"	34 lb	300 lb
3F012230	Flip-Back Desk Arm, Elevating legrest	16"	16"*	18" or 20"	17 ¼"	29"	35½"	24"	12 ¼"	34 lb	300 lb
3F012320	Flip-Back Desk Arm, Swingaway detachable footrest	20"	16"*	18" or 20"	17 ¼"	29"	35½"	28"	12 ¼"	34 lb	300 lb
3F012330	Flip-Back Desk Arm, Elevating legrest	20"	16"*	18" or 20"	17 ¼"	29"	35½"	28"	12 ¼"	34 lb	300 lb
3F012340	Flip-back Full Arm, Swingaway Footrest	20"	16"*	18" or 20"	17 ¼"	29"	35 ½"	28"	12 ¼"	34 lb	300 lb
3F013120	Flip-Back Adjustable Height Desk Arm, Swingaway Footrest	18"	16"*	18" or 20"	17 ¼"	29"	35 ½"	26"	12 ¼"	34 lb	300 lb
3F013220	Flip-Back Adjustable Height Desk Arm, Swingaway Footrest	16"	16"*	18" or 20"	17 ¼"	29"	35 ½"	24"	12 ¼"	34 lb	300 lb
3F013320	Flip-Back Adjustable Height Desk Arm, Swingaway Footrest	20"	16"*	18" or 20"	17 ¼"	29"	35 ½"	28"	12 ¼"	34 lb	300 lb
3F014120	XVI, Flip-Back Desk Arm, Swingaway Footrest	18"	16"*	16.3" or 17.9"	17 ¼"	20"	33"	25.78"	12 ¼"	31 lb	300 lb
3F014220	XVI, Flip-Back Desk Arm, Swingaway Footrest	16"	16"*	16.3" or 17.9"	17 ¼"	20"	33"	23.81"	12 ¼"	31 lb	300 lb
3F014320	XVI, Flip-Back Desk Arm, Swingaway Footrest	20"	16"*	16.3" or 17.9"	17 ¼"	20"	33"	27.75"	12 ¼"	31 lb	300 lb

* + 2" extension

Traveler® L3 Plus – Common Accessories & Parts		HCPCS Code
90763227	Complete Assembly adjustable height desk arm, 1pr	E0973
300-505	Positioning belt 27-47" L (auto-style) push button	E0978
90763062A	Wheel lock lever extension, pair	E0961
90763064A	Heel loop, pair	E0951
90763230A	Anti-tippers adjustable, pair for regular 18" to 20" seat height	E0971
90763233A	Anti-tippers adjustable, pair for 16.3" or 17.9" seat height	
907633CSE	Seat rail plug (dust cap) for seat extension, two pairs per bag	
90763SE16	2" seat extension kit for 16" W	E0981
90763SE18	2" seat extension kit for 18" W	E0981
90763SE20	2" seat extension kit for 20" W	E0981
B10320F	O2-IV holder	E2208
B10325	O2 cylinder holder	E2208
EJO2BAG07	O2 holder bag	E2208