

The Lumex Clinical Care Recliner Wide with Drop Arm has a weight capacity of 450 pounds when evenly distributed. An innovative, easy to actuate Drop-Arm feature allows for a safer and easier lateral transfer for both the patient and caregiver.

Features:

- Dual gas cylinders give smooth power assistance when placing the patient into and out of the Trendelenburg Position
- Wide, 24" seat is more comfortable for the patient
- Maximum Weight Capacity 450 lb – EVENLY DISTRIBUTED
- Two Self Storing Drop Arms are easy to actuate by the caregiver; allows for a safer lateral transfer into the recliner from a bed, wheelchair or stretcher
- Two self storing side tables for patient prep/convenience; smooth surface allows for easy cleaning
- High quality tette 5" swivel casters; three easy to lock total locking casters and one directional caster
- Wall saver armrest tray increases durability of the armrests and protect walls from marring
- Side panels are easily removed to enable access o the recliner for cleaning and maintenance
- Upholstered in healthcare grade vinyl or urethane; available in a broad palette of colors
- Meets California Technical Bulletin CA TB-133 Flammability Test
- All parts can be replaced on recliner greatly enhancing the lifespan of the recliner and reducing its cost in use over time



- Available Options:
 - Heat and Massage Function
 - ◊ Three Heat Settings for patient comfort
 - ◊ Independent Massage Function helps soothe muscles during long treatments
 - Fully Upholstered Armrests
 - Headrest Cover: HRC587 (COLOR CODE)
 - Universal Pillow: FR5659 (COLOR CODE) US

Specifications

Overall Height:49.5" / 125.7 cm
Overall Width:33.0" / 83.8 cm
Overall Depth:44.0" / 111.8 cm
Depth Fully Reclined:71.0" / 180.3 cm
Seat Width:23.5" / 59.7 cm
Seat Depth:22.0" / 55.9 cm
Seat Height:22.0" / 55.9 cm



Upright Position



Upright Position with
Legrest Extended



First Recline Position



Full Recline Position



Trendelenburg Postion