

- Applications: skin or mucous membrane cauterization, removing warts and granulated tissue
- 75% Silver Nitrate and 25% Potassium Nitrate
- Ingredient Net Weight per 6" Applicator tip: 50 mg (37.5 mg Silver Nitrate, 12.5 mg Potassium Nitrate)
- Silver Nitrate tips are soluble in all secretions
- Control strength by dipping in sterile, distilled water
- Rigid wooden applicator for precise placement
- Convenient tube vial
- Made in USA, GF Health Products, Inc. Rhode Island Facility
- ISO 13485 Certified Facility

U.S. Federal Law prohibits dispensing without prescription.

For professional use only.

Distributors (USA) must have an active state wholesale drug license to purchase Grafc® Silver Nitrate Applicators.



Item #	Description	NDC No.	Active Ingredients	UOM
1590	6", Double Dipped	12165-100-01, 12165-100-03	50 mg / applicator	100/tube, 12 tubes/box
1592	12", Double Dipped	12165-200-03	70 mg / applicator	100/tube
1593	18", Double Dipped	12165-300-01	98 mg / applicator	1 ea

