




AUTOMATIC OVERBED TABLES



*Total
Design*
Solution



FINISHES AVAILABLE ON 3DL THERMOFORM OVERBED TABLE TOPS



Aria



Arctic Groovz



Cocobala



Kensington Maple



Shaker Cherry



Solar Oak



Talas Cherry



Wild Cherry



Windsor Mahogany




FINISHES AVAILABLE ON LAMINATE OVERBED TABLE TOPS



Aria




Arctic Groovz




Biltmore Cherry



Cocobala



Empire Mahogany




Fonthill Pear



Kensington Maple




Montana Walnut



Shaker Cherry



Solar Oak



Talas Cherry



Wild Cherry



Windsor Mahogany

AUTOMATIC OVERBED TABLES

Ready Ship - Overbed Table Tops & Base Sets⁵

Item	A793053	A793055	A793020	A793029	A793004	A860053	A860029	A860004	A797053	A797055	A797020	A797029	A797004	A859053	A859029	A859004	A919053	A919055	A919020	A919029	A919004
Low Height ¹	■	■	■	■	■	■	■	■	—	—	—	—	—	—	—	—	—	—	—	—	—
Standard Height ¹	—	—	—	—	—	—	—	—	■	■	■	■	■	■	■	■	—	—	—	—	—
Extended Height ¹	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	■	■	■	■	■
H Base	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	—	■	■	■	■	■
Modified H Base	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	—	—	—	—	—
3DL ²	—	—	—	—	—	■	■	■	—	—	—	—	—	■	■	■	—	—	—	—	—
Laminate ³	■	■	■	■	■	—	—	—	■	■	■	■	■	—	—	—	■	■	■	■	■
Flat Edge	■	■	■	■	■	—	—	—	■	■	■	■	■	—	—	—	■	■	■	■	■
Raised Edge	—	—	—	—	—	■	■	■	—	—	—	—	—	■	■	■	—	—	—	—	—
Solar Oak	■	—	—	—	—	■	—	—	■	—	—	—	—	■	—	—	■	—	—	—	—
Talas Cherry	—	■	—	—	—	—	—	—	—	■	—	—	—	—	—	—	—	■	—	—	—
Cocobala	—	—	■	—	—	—	—	—	—	—	■	—	—	—	—	—	—	—	■	—	—
Windsor Mahogany	—	—	—	■	—	—	■	—	—	—	—	■	—	—	■	—	—	—	—	■	—
Wild Cherry	—	—	—	—	■	—	—	■	—	—	—	—	■	—	—	■	—	—	—	—	■
Weight Capacity EVENLY DISTRIBUTED	lb	40	40	40	40	40	40	40	40	40	40	40	40	40	40	40	100	100	100	100	100
UOM		1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea	1/ea

■ Standard Feature ■ Optional Feature ■ Available Style ■ Available Finish — Not Available/Does Not Apply

Made To Order
Top Options:⁷

Item	3DL ² Laminate ³ material	Width inch	Length inch	Raised Edge Flat Edge top profile	Finish	UOM
A28B Laminate Top	— ■	15	30	— ■	Specify ^{3,4}	1/ea
A28-3DL 3DL OBT Top	■ —	15	30	— ■	Specify ^{2,4}	1/ea
A28-REB 3DL Raised Edge	■ —	15	30	■ —	Specify ^{2,4}	1/ea
A28A Lam. Oversized	— ■	18	35	— ■	Specify ^{3,4}	1/ea

Base
Options:⁷

Item	Height	Width inch	Length inch	Height Range ¹ inch - low inch - hi	Weight Capacity EVENLY DISTRIBUTED lb	UOM
A793000	Low	15.30	28.50	21.00 32.00	40	1/ea
A797000	Standard	15.30	28.50	27.75 40.00	40	1/ea
A919000	Extended	18.20	29.30	29.90 43.00	100	1/ea

Accessories:⁶

Item	UOM
A28VAN Vanity with Pop Up Mirror [FITS ALL]	1/ea

¹ height adjusts infinitely
² 3DL finish swatches on page 2
³ laminate finish swatches on page 3

⁴ must specify OBT top finish when ordering
⁵ Ready Ship tops & bases sold together
⁶ available for either Ready Ship or Made to Order
⁷ Made to Order tops and bases each sold separately

• additional styles of OBT tops may be available upon request
• all measurements include a manufacturing variance of ± 0.25 inch
• cover images - left: A859004 40 lb OBT right: A919055 100 lb OBT



WARNINGS: Cancer and Reproductive Harm -
www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company



MDSAP
MEDICAL DEVICE SUPPLIER ASSESSMENT PROGRAM

© 2020, GF Health Products, Inc. All Rights Reserved.
Total Design Solutions is a service mark of GF Health Products, Inc. Graham-Field
and Basic American Medical Products are trademarks of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website
www.grahamfield.com

**graham
field**

GF120004RevG20