



BLADES

STERILE DISPOSABLE MINI SCALPELS

Feather sterile disposable mini scalpels feature a slender plastic handle that is ideal for surgery in limited areas such as the brain, ear, nose, mouth, blood vessels, or nerves, as well as in plastic or neonatal surgery.

Feather, a blade manufacturer of global renown, utilizes high-precision technology and carefully selected materials. Each Feather scalpel is made in Japan and meets strict ISO standards, resulting in exceptionally high quality and consistent strength and sharpness.

Features

- Sterile: Gamma-sterilized and ready for use immediately after opening
- Single-Use/Disposable: Helps prevent patient cross-contamination
- Heat-treated stainless steel blade: Double-honed blade ensures sharp cutting edges with superb durability
- Individually packaged: Hermetically sealed aluminum foil packaging complete with clear window for product selection verification and pull tabs for easy opening
- ABS* handle includes textured grip
- Protective snap-off cover: Protects blade and user until ready to use
- Available in 3 different sizes
- Color-coded by size: Packaging assists with blade identification
- Individually packaged 20 per dispensing box



Item No.	Sterile, Single-Use, Stainless Steel Blade with (ABS*) Plastic Handle	Blade Angle	Color Code	UOM
2978#11		11		20/box dispenser
2978#14		14		20/box dispenser
2978#15		15		20/box dispenser

Color variations may occur due to monitor or printer settings.



GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company
MDSAP
MEDICAL DEVICE SUPPLY PROGRAM

© 2021, GF Health Products, Inc. All Rights Reserved.
GF Health Products, Inc. is the US distributor for FEATHER
FEATHER is a registered trademark of FEATHER Safety Razor Co., Ltd



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com



NOTE: Dispose of scalpel in an approved sharps container in accordance with your facility's sharps disposal procedures

Product Hyperlinks Enabled

* Acrylonitrile Butadiene Styrene

graham
field

GF2100013 Rev/B21