



BLADES

SAFESHIELD STERILE DISPOSABLE SCALPELS

Designed with an ergonomically rounded SUREgrip™ plastic handle that fits the hand comfortably, and outfitted with a stainless steel blade, Safeshield™ Safety Scalpels feature a locking retractable protective shield that helps prevent injuries before, during and after use.

Feather, a blade manufacturer of global renown, utilizes high-precision technology and carefully selected materials. Each Feather scalpel is made in Japan and meets strict ISO standards, resulting in exceptionally high quality and consistent strength and sharpness.

Features

- Single Use/Disposable: Helps prevent patient cross-contamination
- Individually packaged: Non-woven Tyvek cloth package complete with clear window for product selection verification and pull tabs for easy opening
- Clear retractable safety shield: Can be operated with one hand
- Auditory confirmation: Transparent blade shield makes audible click when locked in position with blade covered or uncovered
- Available in 7 sizes: Including most common blade sizes and select specialty blade sizes
- Color-coded by size: Release tab on protective shield and coordinating packaging assist with blade identification
- Ruler printed on handle: Includes cm and mm markings
- Similar weight as conventional disposable scalpel
- Sterile: Gamma-sterilized and ready for use immediately after opening
- Individually packaged 10 per dispensing box

Item No.	Stainless Steel Blades	Blade Size	Color Code	UOM
2980#10		10	Orange	10/box dispenser
2980#11		11	Pink	10/box dispenser
2980#14		14	Red	10/box dispenser
2980#15		15	Purple	10/box dispenser
2980#20		20	Gold	10/box dispenser
2980#21		21	Green	10/box dispenser
2980#22		22	Black	10/box dispenser



Reasons to Make the Switch to Safety Scalpels

- Industry requirements and professional organizations demand them
- Protect employees from on-the-job injuries
- Create a safe place to work; have a vested interest in avoiding the health, staffing, and legal costs of sharp injuries

Color variations may occur due to monitor or printer settings.



GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company

© 2021, GF Health Products, Inc. All Rights Reserved.
GF Health Products, Inc. is the US distributor for FEATHER
FEATHER is a registered trademark of FEATHER Safety Razor Co., Ltd



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com



NOTE: Dispose of scalpel in an approved sharps container in accordance with your facility's sharps disposal procedures

Product Hyperlinks Enabled

graham
field

GF-1500081 Rev C21