

## THERMOMETERS

# DIGITAL THERMOMETER, °F/°C, FLEXIBLE TIP, 24 BULK PACKAGE

### Features

- Auto Shut Off
- Easy to read LCD Display
- Records and Displays oral, rectal, axillary temperature



### Specifications:

	Width		Length		Height		Range		Accuracy		Memory	Battery Life	Product Weight		UOM
	inch	cm	inch	cm	inch	cm	F°	C°	F°	C°			oz	g	
1858	0.7	1.8	4.9	12.4	0.35	0.9	90 - 109.9	32.22 - 43.28	+/- 0.2°, 95.9°-107.6°	+/- 0.11°, 35.5°-42°	Last measured Value	Approximately 200 hours	~.35 including battery	~ 10 including battery	24/bx
1858CASE									+/- 0.4° under 95.9° or over 107.6°	+/- 0.13° under 35.5° or over 42°					240/cs

**WARNINGS:** Cancer and Reproductive Harm - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Color variations may occur due to monitor or printer settings.



Intertek

GF Health Products, Inc. is an  
ISO 13485:2016  
MDSAP Certified Company



© 2020, GF Health Products, Inc. All Rights Reserved.  
Gendron is a registered trademark of GF Health Products, Inc.



+1770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)

**graham field**

GF1000004RevC21