


REPLACEMENT TOPS



REPLACEMENT TOPS

Make your exam tables look new at a fraction of the cost. All components within in the replacement top are all new including wood, foam, vinyl and hardware.

- Industry standard colors *plus*
- Replacement Tops fit leading brand tables



Standard Color Options (-xx)

Available in 24 color choices to complement most decors.

WH Super White	DO Doe	PM Palm	SH Shadow	AR Artichoke	TE Teal
NO Northwoods	DB Dusty Blue	PU Plum	WW Wedgewood	BT British Tan	GR Gray
AL Alabaster	ST Saddle Tan	PL Pearl	TL Tea Leaf	MO Moss	OG Ocean Gray
AQ Aqua Green	BU Blue	VI Violet	NA Navy	WG Warm Gray	BL Black

Switch Card is available from your representative

WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP) which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to - www.P65Warnings.ca.gov/furniture

Color variations may occur due to monitor or printer settings.



© 2021, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com

Specifications:

	100 104		204 ¹ 223 ²		304	
	inch	cm	inch	cm	inch	cm
Overall Width:	27.50	69.85	27.50	69.85	27.50	69.85
Overall Length:	55.00	139.70	55.00	139.70	55.00	139.70
Footrest Width:	20.00	50.80	19.00	48.26	20.00	50.80
Footrest Length:	18.00	45.72	18.00	45.72	18.00	45.72
	lb	kg	lb	kg	lb	kg
Shipping Weight:	55	24.95	55	24.95	55	24.95

Features:

Velcro® Attachment:	—	—	—
Paper Roll Holder:	included	—	—
Mounting Hardware:	included	—	—
Lift Handles:	—	not included ³	—
Top Replacement Knobs:	—	not included ³	—
Individual Product Links:	link	link	link

¹ If OEM model is seamless top, request paper roll holder when ordering

² For OEM models manufactured after 2009 only

³ Lift Handles and Top Attachment Knobs for exam table year models 2009 and after are not included

Specifications:

	404 405		75L 411		420 460	
	inch	cm	inch	cm	inch	cm
Overall Width:	27.50	69.85	27.50	69.85	27.50	69.85
Overall Length:	60.00	152.40	70.00	177.80	60.00	152.40
Footrest Width:	20.50	52.07	19.00	48.26	19.00	48.26
Footrest Length:	17.50	44.45	18.00	45.72	18.00	45.72
	lb	kg	lb	kg	lb	kg
Shipping Weight:	35	15.88	35	15.88	55	24.95

Features:

Velcro® Attachment:	included	included	—
Paper Roll Holder:	—	—	included
Mounting Hardware:	—	—	included
Lift Handles:	—	—	—
Top Replacement Knobs:	—	—	—
Individual Product Links:	link	link	link

— Denotes Feature Not Available / Does Not Apply

Product Hyperlinks Enabled