

EXAM TABLES



EXAM TABLES


The medical exam table is essential in any health care facility. Hausted provides high-quality options to meet the needs of hospitals, clinics, ERs and more.

- Five Position Continuous Lock Stirrups, the most positions in the industry
- Pull-Out Leg Rest with Integrated Treatment Tray
 - ♦ Tray is independent of the Legrest and may be accessed easily at any time
- An all-metal construction for the base table provides longevity and recyclability
- Powder Coat Finish - Foghorn
- Pelvic Tilt



4200 Series

4200-xx exam table with three pass-through drawers, two storage drawers and stirrups



4600 Series

4600-xx high/low power exam table with power assisted back, hand control and stirrups. UL listed (60601-1) DC electromechanical drive. Hand pendant controls hi-low, back, flat, return to chair



Standard Color Options (-xx)


Available in 24 color choices to complement most decors.

WH Super White	DO Doe	PM Palm	SH Shadow	AR Artichoke	TE Teal
NO Northwoods	DB Dusty Blue	PU Plum	WW Wedgewood	BT British Tan	GR Gray
AL Alabaster	ST Saddle Tan	PL Pearl	TL Tea Leaf	MO Moss	OG Ocean Gray
AQ Aqua Green	BU Blue	VI Violet	NA Navy	WG Warm Gray	BL Black

Switch Card is available from your representative

WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP) which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to - www.P65Warnings.ca.gov/furniture

Color variations may occur due to monitor or printer settings.



© 2021, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com

Specifications:

	4200 SERIES		4600 SERIES	
	inch	cm	inch	cm
Overall Width:	27.00	68.58	27.00	68.58
Overall Length: <i>legrest extended</i>	72.00	182.88	72.00	182.88
Overall Height: <i>minimum</i>	32.00	81.28	18.50 ± 0.5	46.99 ± 1.27
Overall Height: <i>maximum</i>	32.00	81.28	35.50 ± 0.5	90.17 ± 1.27
Foot Step Width:	18.00	45.72	—	—
Foot Step Length:	17.00	43.18	—	—
Legrest Width:	18.00	45.72	18.00	45.72
Legrest Length:	17.00	43.18	17.00	43.18
Paper Roll Holder	18.00	45.72	18.00	45.72
Widths Available:	19.00	48.26	19.00	48.26
	20.00	50.80	20.00	50.80
Drawer Width: ¹	18.00	45.72	—	—
Drawer Length: ¹	16.50	41.91	—	—
Drawer Height: ¹	4.00	10.16	—	—
Foot Step Capacity: <i>EVENLY DISTRIBUTED</i>	800	362.87	—	—
Weight Capacity: <i>EVENLY DISTRIBUTED</i>	500	226.80	700	317.51
Shipping Weight:	335	151.95	395	179.17

Features:

Base:	pass-through drawer	dual pedestal
Stirrups:	five position continuous lock	five position continuous lock
Treatment Tray:	integrated	integrated
Powder Coat Finish:	foghorn	foghorn
Pull Out Leg Rest:	standard	standard
Pelvic Tilt:	standard	standard

Options:

	available at additional cost	
	order number	order number
None:	4200-xx	4600-xx
Drawer Warmer: <i>w/ integrated electrical outlet</i>	4200-DW-xx	—
Drawer Warmer* <i>includes physician storage</i>	—	4600-PS-DW-xx
Integrated Outlet:	4200-E0-xx	—
Physician Storage* <i>without drawer warmer</i>	—	4600-PS-xx
Assist Arms*	—	4600-AA-xx
Foot Switch*	—	4600-FS-xx
Fully Loaded: <i>options marked *</i>	—	4600-FL-xx
Replacement Top:	GF420RT-xx	GF460RT-xx

* Options Included in the Fully Loaded Package

¹ Interior Dimensions

Product Hyperlinks Enabled