



EXTENDED CARE RECLINER

PREFERRED CARE[®] RECLINER

Infinite recliner back patented mechanism provides infinite positioning, ranging from 95 degrees to 145 degrees for increased resident comfort and pressure relief. The reclining back can be locked in any position by the caregiver.

Reduced armrest and seat height enables the recliner to fit under dining/activity tables and allows resident to "scoot" and more easily reposition chair by themselves.

Extra-wide padded armrests and side panels maximize comfort.

Two front swivel casters and two rear locking casters enable the recliner to be easily maneuvered in small spaces and moved down hallways.

A unique suspension system self-adjusts to each person, cradling the resident from head to heel and virtually eliminating the risk of "bottoming out".

Quick-release upholstery allows for easy cleaning, removal and replacements. Vinyl upholstered back is sulfide stain, mildew and oil resistant.

Optional Accessories:

- Anti-Tippers
- IV Pole & Mount
- Headrest Pillow
- Head Support
- Body Bolsters
- Adjustable Footrest
- Foot-Drop Assembly
- Activity Tray Table



Blue Ridge (427)

Jade (857)

Rosewood (863)

Model (xxxx represents fabric color)	Standard Width Seat	Extra-Wide Seat	Drop-Arm	Tilt-in-Space	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from Seat		Width between arms		Depth fully reclined		Weight Capacity EVENLY DISTRIBUTED		UOM
					inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
					FR565Gxxxx	●				28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	20	50.8	8	20.32	19	
FR565DGxxxx	●	●			28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	22	55.88	8	20.32	19.5	49.53	66	167.64	250	113.4	1 ea
FR565TGxxxx	●		●		28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	22	55.88	8	20.32	19.5	49.53	66	167.64	250	113.4	1 ea
FR565WGxxxx		●			32	81.28	37	93.98	48	121.92	24	60.96	19.5	49.53	20	50.8	8.5	21.59	24	60.96	66.5	168.91	400	181.44	1 ea

Color variations may occur due to monitor or printer settings.