



CANES

STANDARD WOODEN CANES

Features

- A classic style made of select hardwood that is satin-smooth to the touch
- Available in a variety of attractive finishes
- Comes complete with rubber tip
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		Diameter		Length		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5180A	Walnut Finish	0.875	2.22	36	91.44	250	113.40	6/cs
5181A	Mahogany Finish	0.875	2.22	36	91.44	250	113.40	6/cs
5182A	Walnut Finish	1	2.54	36	91.44	250	113.40	6/cs
5183A	Mahogany Finish	1	2.54	36	91.44	250	113.40	6/cs
5184A	Walnut Finish, Extra Long	1	2.54	42	106.68	250	113.40	6/cs



Intertek

GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company



© 2015, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.



+1.770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Color variations may occur due to monitor or printer settings.



GF-15-0005-RevB21