



CANES

STYLUS
OFFSET CANES

Features

- Available in a variety of attractive finishes
- Manufactured with sturdy, extruded aluminum tubing
- Offset handle places the user's weight directly over the base for maximum balance and control
- Soft, contoured foam grip for maximum comfort
- Easy push-button height adjustment from 30" - 39" with locking collet for added security
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5950SB	Blue	30	76.2	39	99.06	250	113.4	6/cs
5950SBK	Black	30	76.2	39	99.06	250	113.4	6/cs
5950SL	Lavender	30	76.2	39	99.06	250	113.4	6/cs
5950SR	Red	30	76.2	39	99.06	250	113.4	6/cs



GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company



© 2012, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.



+1.770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Color variations may occur due to monitor or printer settings.

