

## CLINICAL CARE RECLINERS

# CLINICAL CARE RECLINER - WIDE

Three positions for patient comfort. The fourth position is a caregiver-activated Trendelenburg (or Shock) Position.

Improved wider and softer armrest for patient comfort.

Quick Ship available.

Dual gas cylinders give smooth power assistance during Trendelenburg activation.

Wide, ergonomically designed padded armrests for improved patient comfort.

Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with one directional locking caster.

Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Side table surface is smooth and features tapered edges for easy cleaning.

Wide selection of fabrics and colors. See the [Swatch Card](#) for details.

### Options:

- Heat and Massage
- Heat Only feature
- Fully Upholstered Armrests
- IV Pole Holder / IV Pole
- Headrest Cover
- Universal Pillow
- Drop-Arm (FR587WD)



Multiple Positions for Enhanced Patient Comfort

Model (xxxx represents fabric colors)	Wall Saver Armrest Edges	Fully Upholstered Arms	Heat & Massage	Heat Only	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from Seat		Width between arms		Width between side tables in upright position		Depth fully reclined (Trendelenburg)		Weight Capacity EVENLY DISTRIBUTED		UOM
					inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
FR587W	●				33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WH	●		●		33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WHO	●			●	33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WHxxxxPA	●				33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WHxxxxPA	●	●			33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea
FR587WHOxxxxPA	●		●		33	83.8	44	111.8	49.5	125.7	23.5	59.7	20.5	52.1	22	55.9	7	17.8	23	58.4	53	134.6	80	203.2	450	204.1	1 ea

Color variations may occur due to monitor or printer settings.