



NON- FOLDING WALKERS

FREEDOM WALKER

Features

- Encourages independence and self esteem by allowing freedom of movement and exercise
- Restraint free opening/closing gate with a safety knob
- Mesh available in 13 colors
- Adjustable safety belt to help prevent injuries
- Solid base seat with white cushion for extra comfort
- Constructed of the highest quality health care grade PVC
- Maximum Weight Capacity
300 lb - EVENLY DISTRIBUTED

Royal Blue mesh will be provided for online orders.

If you desire another color, please order via Customer Service at 1-800-347-5678.



80500	internal Width		external Width		minimum Seat Height		maximum Seat Height		minimum Arm Height		maximum Arm Height		overall Height		overall Length		Walking Space front of seat to gate		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
	18	45.72	30.5	77.47	20.25	51.44	24.75	62.87	28.75	73.03	33.25	84.46	40.75	103.51	35	88.9	15	38.10	300	136.08	2/cs