



## FOLDING WALKERS

# DUAL RELEASE FOLDING WALKERS

### Features

- Durable 1" aluminum tubing provides maximum strength while remaining lightweight
- Raised H-frame design facilitates over-toilet positioning
- Contoured, textured handgrips offer a safe, secure hand hold
- Dual release folding mechanism enables user to fold walker sides independently
- Adjustable height fits a broad range of user heights
- Accommodates all leg attachments
- Folded depth 4"
- Designed to meet the requirements of HCPCS code: (without wheels) E0135 (with wheels) E0143
- Includes laminated color operating instructions
- Models with wheels come with rear glide caps to facilitate use on all surfaces
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		5" Fixed Wheels	Width		Width		overall		Depth		Depth		minimum		maximum		user		Weight		UOM
			inside handgrips	inside rear legs	Width	Depth	folded	Height	Height	Height Range	Capacity										
		No wheels	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg			
616270A-4	Adult	●	16.5	41.91	20.125	51.12	26	66.04	16.5	41.91	4	10.16	34.5	87.63	41.5	105.41	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs
616370A-4	Junior	●	16.5	41.91	18.625	47.31	25.38	64.45	15.5	39.37	4	10.16	29.5	74.93	34.5	87.63	4' 7" - 5' 3"	1.4 - 1.6	300	136.08	4/cs
616070A	Adult	●	16.5	41.91	20.125	51.12	22.75	57.79	16.5	41.91	4	10.16	32.5	82.55	41.5	105.41	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs



Intertek

GF Health Products, Inc. is an  
ISO 13485:2016  
MDSAP Certified Company



© 2015, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)

Color variations may occur due to monitor or printer settings.

graham  
field

GF-15-00016RevC21