



FOLDING WALKERS

TWO BUTTON FOLDING WALKER  
WITH 5" WHEELS

**Features**

- Easy push dual-release folding mechanism enables user to fold walker sides independently
- Sturdy 1" diameter steel construction ensures maximum strength
- PVC contoured hand grip
- Comes with 5" wheels (Replacement leg set with wheels item #716280-LEGS)
- Rear glide tips allow for walker to move easily and smoothly over most surfaces
- Designed to meet the requirements of HCPCS code: E0143
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- One (1) year limited warranty

**Replacement Parts**

		UOM
716280-LEGS	Replacement Leg Set	1 set
716280-TIPS	Replacement Tips	1 pr



		Width inside handgrips		overall Width		Depth front to rear		minimum Height		maximum Height		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
716280A-4	Two Button Folding Walker	17.5	44.45	24.5	62.23	20.25	51.44	32.5	82.55	39.5	100.33	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs