

DIAGNOSTIC INSTRUMENTS

NEUROLOGICAL SET

A neurological examination assesses motor and sensory skills, hearing and speech, vision, coordination, and balance. It may also test mental status, mood, and behavior. The examination uses tools such as reflex hammers and a pinwheel.

Each neurological set includes:

7" Buck Neurological Hammer

- Octagonal chrome plated handle
- Dual latex-free heads
- Removable needle and brush

7.5" Taylor Hammer

- Latex-free black triangular head
- Chrome plated handle

8" Babinski Percussion Hammer

- Adjustable latex-free head can be oriented in 2 directions
- Octagonal chrome plated handle

7" Wartenburg Pinwheel

- Chrome plated handle and rotating spur

Also includes

- Additional babinski head
- Additional 12" handle for babinski hammer

Set comes complete in a blue 13" x 1" (33 x 2.5 cm) container

Latex Free



		UOM
1303	Neurological Set	1 set



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.

graham
field