

EXAM ROOM SUPPLIES

LABELED GLASS SUNDRY JARS

Grafco[®] Labeled Glass Sundry Jars are excellent for storing and dispensing a variety of supplies, making them ideal for hospitals and clinics.

Features

- Transparent flint glass construction provides easy visual of contents inside
- Each jar permanently imprinted with blue labels: Gauze, Applicators, Tongue Depressors, Cotton and Bandages
- Reusable
- Not autoclavable
- Available in sets including four (4) or five (5) flint glass jars complete with removable overlapping polished aluminum lids
- Latex free
- One (1) year limited warranty



		Height		Diameter		case Weight		UOM
		inch	cm	inch	cm	lb	kg	
3454	Set of 5	7	17.8	4.25	10.8	7.5	3.4	5 ea/st
3454A A	Applicators Jar w/ lid	7	17.8	4.25	10.8	1.44	0.65	1 ea
3454A B	Bandages Jar w/ lid	7	17.8	4.25	10.8	1.44	0.65	1 ea
3454A C	Cotton Jar w/ lid	7	17.8	4.25	10.8	1.44	0.65	1 ea
3454A G	Gauze Jar w/ lid	7	17.8	4.25	10.8	1.44	0.65	1 ea
3454A T	Tongue Depressors w/ lid	7	17.8	4.25	10.8	1.44	0.65	1 ea
3454-4	Set of 4 (excluding tongue depressors jar)	7	17.8	4.25	10.8	7.125	3.23	4 ea/st
3454C	Replacement Stainless Steel Lid	7	17.8	4.25	10.8	—	—	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2016, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.

graham
field

GF1600114Rev/D22