





EXAM ROOM SUPPLIES

SILVER NITRATE APPLICATORS

- Applications: skin or mucous membrane cauterization, removing warts and granulated tissue
 75% Silver Nitrate and 25% Potassium Nitrate
 Ingredient Net Weight per 6" Applicator tip: 50 mg (37.5 mg Silver Nitrate, 12.5 mg Potassium Nitrate)

- Silver Nitrate tips are soluble in all secretions
- Control strength by dipping in sterile, distilled water
 Rigid wooden applicator for precise placement
 Convenient tube vial

- Made in USA, GF Health Products, Inc. Rhode Island Facility
- ISO 13485 Certified Facility









For professional use only. Distributors (USA) must have an active state wholesale drug license to purchase Grafco® Silver Nitrate Applicators.

		Size		Active Ingredients	Double Dipped	
Item	NDC N°	inch	cm	mg	2	UOM
1590	12165-100-01	6.00	15.24	50		100/tube, 12 tubes/box
1592	12165-200-03	12.00	30.48	70		100/tube
1593	12165-300-01	18.00	45.72	98		1/ea
	■ Denotes Standard Feature				— Not Available / Does Not Apply	



«GF Health Products, Inc. is an intertek ISO 13485:2016 / MDSAP Certified Company











