

## DAILY LIVING AIDS

# DELUXE SOCK AID

Essential item to help avoid bending or twisting after hip or knee surgery

### Features

- Textured cloth side prevents sock from slipping
- Smooth nylon side reduces friction with foot
- Large easy to grasp loop handles
- 29" long handle straps are ideal length
- Attractive Retail Box



	overall Length		UOM
	inch	cm	
5592R-1	29	73.66	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF2000183RevB22