


DAILY LIVING AIDS

GAIT BELT

Essential item to help lift and transfer a patient and to help reduce caregiver injury.

Features

- Available in 54" or 72" lengths
- Strong and durable 2" webbing
- Durable metal buckle with teeth to better secure webbing
- Attractive Retail Box



	overall Length		UOM
	inch	cm	
5597-54R-1	54	137.16	1 ea
5597-72R-1	72	182.88	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF2000185RevB22