

## ROLLATORS

## HYBRIDLX ROLLATOR TRANSPORT CHAIR

The HybridLX is an innovative rollator that combines the function of a wheeled walker and a transport chair into one lightweight design. Keeping the patient's desire for independence and comfort in mind, the HybridLX is designed with secure handgrips, a comfortable and wide seat and backrest, and the effortless design to transition from a rollator to a transport chair.

## **Features**

- Lightweight frame durable, aluminum frame weighs only 22 lbs (this includes footrests) making it easier for users to move about freely
- Hybrid concept designed to be used primarily as a rollator while easily transforming into a transport chair with removable footrests (adjustable in length) when user needs a short rest
- Ultra-Compact design easy to assemble without the need for tools and folds effortlessly for convenient transport and storage
- Easy-to-use handbrakes allow for short and longer period braking when used as rollator and transport chair (user must be accompanied when using as a transport chair)
- Flexible, soft polyester backrest and 18" wide mesh seat with pull up strap for easy folding provides comfort for user
- Large 8" wheels for easy rolling indoor and outdoor on rough surfaces
- Complete with storage pouch for personal items
- Available in attractive, contemporary colors: Titanium and Majestic Blue
- Designed to meet the requirements of HCPCS Codes: E0143 + E0156
- Three (3) year limited warranty



	overall <b>Width</b>			overall <b>Depth</b> with footrests		Depth front to rear			Width inside handgrips			minimum Handle Height		<sub>maximum</sub> Handle Height			Wheel diameter		
	inch	cm		inch	cm		inch	cm		inch	cm		inch	cm	inch	cm	17	inch	cm
LX1000T Titanium	26	66.04		40	101.60		28.5	72.39		18.5	46.99		30	76.20	37	93.98		8	20.32
LX1000B Majestic Blue	26	66.04		40	101.60		28.5	72.39		18.5	46.99		30	76.20	37	93.98		8	20.32

	Seat '	Width	Seat	Depth	Seat I	Height	Height	Product Weight			Ca EVENLY			
	inch	cm	inch	cm	inch	cm	feet	m		lb K	g	lb	Kg	UOM
LX1000T Titanium	18	45.72	11	27.94	21.25	53.98	5' 4" - 6' 2"	1.63 - 1.88	В	22 9.9	98	300	136.08	1 ea
LX1000B Majestic Blue	18	45.72	11	27.94	21.25	53.98	5' 4" - 6' 2"	1.63 - 1.88	8	22 9.9	98	300	136.08	1 ea





GF Health Products, Inc. is an intertek ISO 13485:2016 / MDSAP Certified Company











