

## FOLDING WALKERS

# WALKANE

### Features

- Lighter than a walker; more stable than a cane
- Folds easily with one hand
- Ideal for users who do not require a walker but need a wider base for support and for individuals with the use of only one hand or arm
- Contoured, foam hand grips are comfortable and provide secure hand hold
- Comfortable, dual level handgrips provide support from a seated to standing position
- Folds to compact size for easy storage and transport
- Comfortable, dual level handgrips provide support from a seated to standing position
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0135
- Maximum Weight Capacity  
300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	Width inside rear legs		overall Width		Depth front to rear		minimum Height		maximum Height		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
6015A	13.25	33.66	16.5	41.91	13.25	33.66	29.5	74.93	34.5	87.63	5' 1" - 5' 10"	1.55 - 1.78	300	136.08	2/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

graham  
field

GF0900240 Rev E22