

## SHOWER TRANSPORT CHAIRS


# PVC SHOWER CHAIR/COMMODE

### Features

- Constructed of the highest quality health care grade PVC
- Deluxe elongated extra heavy-duty, contoured open from seat
- Dual reinforcement into the caster fittings adds additional support
- Quick drying, removable mesh back support is standard in all shower chairs
- Available in 3 anti-bacterial mesh or anti-bacterial solid vinyl colors
- Royal Blue mesh is the default color shipped. If you desire another color or vinyl, please order via Customer Service at 770-368-4700
- Ergonomic curved arms provide support without sharp edges and include anti-slip handgrips for added safety
- Threaded stem casters for additional safety
- Optional telescoping footrest feature allows footrest to move in or out with ease
- Optional positioning belt available separately. (Item#89301)
- One (1) year limited warranty on PVC frame; six (6) month limited warranty on all other materials



### Standard Mesh and Solid Vinyl Colors



		overall Height		overall Width		overall Depth		Seat Height		Seat Width		Seat Depth		Width between arms		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
89110	with foot rest, 3" twin wheel	40	101.6	22	55.88	21	53.34	21	53.34	18	45.72	18	45.72	18	45.72	300	136.08	1 ea
89150	without foot rest, 3" twin wheel	40 - 45	101.6 - 114.3	22	55.88	21	53.34	20-25	50.8 - 63.5	18	45.72	18	45.72	18	45.72	250	113.4	1 ea
89210	with foot rest, 3" heavy duty caster	40	101.6	26	66.04	23	58.42	21	53.34	22	55.88	19.5	49.53	22	55.88	375	170.1	1 ea
89250	without foot rest, 4" heavy duty caster	40	101.6	30	76.2	30	76.2	20	50.8	26	66.04	22	55.88	26	66.04	425	192.78	1 ea
89251	with foot rest, 4" heavy duty caster	40	101.6	30	76.2	30	76.2	20	50.8	26	66.04	22	55.88	26	66.04	425	192.78	1 ea
89350*	without foot rest, 5" heavy duty caster	43	109.22	34	86.36	28.25	71.76	23	58.42	30	76.20	21.5	54.61	30	76.20	700	317.51	1 ea
89351*	with foot rest, 5" heavy duty caster	43	109.22	34	86.36	28.25	71.76	23	58.42	30	76.20	21.5	54.61	30	76.20	700	317.51	1 ea

\*89350 and 89351 cannot be used over a toilet



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF0700131RevD22