



BARIATRIC HOMECARE BED

HOME CARE 3648 SERIES

**Features**

- Two-piece frame construction for ease of delivery and transport
- Solid pan-style sleep deck
- Deluxe caster set features durable Total-Lock and Steer-Lock casters
- Fully electric, pendant control
- Standard 3-Function Positioning hand control
- Head end side rails (shown on bed at right) are included
- 650 lb (194.8 kg<sup>1</sup>) weight capacity<sup>1</sup> –EVENLY DISTRIBUTED

**3648** Home Care Bed Product Info



| Options:             | Home Care Bed Only                  | Split Frame                         | Tubular Side Rails                  | Head & Foot Boards                  |
|----------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| GS3648               | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> |
| GS3648B <sup>2</sup> | <input checked="" type="checkbox"/> | <input type="checkbox"/>            | <input type="checkbox"/>            | <input type="checkbox"/>            |

**Specifications:**

| Series | Bed Deck Width |        | Bed Deck Length |        | Bed Deck Height Range |                 | Weight <sup>1</sup> Capacity |                 |
|--------|----------------|--------|-----------------|--------|-----------------------|-----------------|------------------------------|-----------------|
|        | inch           | cm     | inch            | cm     | minimum<br>inch       | maximum<br>inch | lb                           | kg <sup>1</sup> |
| 3648   | 48.00          | 121.92 | 84.00           | 213.36 | 15.00                 | 38.10           | 650                          | 294.84          |

\* options available at additional cost † standard rounding method used <sup>1</sup> maximum weight capacity - evenly distributed - refer to user manual <sup>2</sup> accessories sold separately



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.