



Vista® IC

The industry compatible Vista IC incorporates the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction and the latest industry standard front rigging, upholstery mount spacing, and detachable armrest (with select models) to provide the most effective combination K0001/K0002 wheelchair available.

Standard Features

- Industry Compatible front rigging that interchanges with other brands.
- Industry Compatible arm sockets accept the latest version of other popular models.
- Industry Compatible hole spacing for seat and back upholstery allows use of many after-market products.
- Two position axle caster fork provides full 2" seat height change.
- Expanded applications with height adjustable arms and seat rail extension kits.
- 250 lb weight capacity for 16" and 18" models.
- 300 lb weight capacity for 20" models.

Warranty

- Limited Lifetime on side frames and crossbars.



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com



Footrest interchangeable with other brands of wheelchairs



Dual axle position in lower position.

Vista® IC (Industry Compatible)		Dimensions and Weights								
		Seat			Back	Overall			Wt.*	
53010100	Fixed Full Arm, Swingaway Footrest	18"	16"	19 3/4"	16 1/2"	41"	36 1/2"	27 1/2"	12"	39 lb
53010110	Fixed Full Arm, Elevating Legrest	18"	16"	19 3/4"	16 1/2"	41"	36 1/2"	27 1/2"	12"	39 lb
53010120	Detachable Desk Arm, Swingaway Footrest	18"	16"	19 3/4"	16 1/2"	41"	36 1/2"	27 1/2"	12"	39 lb
53010130	Detachable Desk Arm, Elevating Legrest	18"	16"	19 3/4"	16 1/2"	41"	36 1/2"	27 1/2"	12"	39 lb
53010140	Detachable Full Arm, Swingaway Footrest	18"	16"	19 3/4"	16 1/2"	41"	36 1/2"	27 1/2"	12"	39 lb
53010150	Detachable Full Arm, Elevating Legrest	18"	16"	19 3/4"	16 1/2"	41"	36 1/2"	27 1/2"	12"	39 lb
53010160	Adjustable Height Desk Arm, Swingaway Footrest	18"	16"	19 3/4"	16 1/2"	41"	36 1/2"	27 1/2"	12"	39 lb
53010170	Adjustable Height Desk Arm, Elevating Legrest	18"	16"	19 3/4"	16 1/2"	41"	36 1/2"	27 1/2"	12"	39 lb
53010200	Fixed Full Arm, Swingaway Footrest	16"	16"	19 3/4"	16 1/2"	41"	36 1/2"	25 1/2"	12"	39 lb
53010210	Fixed Full Arm, Elevating Legrest	16"	16"	19 3/4"	16 1/2"	41"	36 1/2"	25 1/2"	12"	39 lb
53010220	Detachable Desk Arm, Swingaway Footrest	16"	16"	19 3/4"	16 1/2"	41"	36 1/2"	25 1/2"	12"	39 lb
53010230	Detachable Desk Arm, Elevating Legrest	16"	16"	19 3/4"	16 1/2"	41"	36 1/2"	25 1/2"	12"	39 lb
53010240	Detachable Full Arm, Swingaway Footrest	16"	16"	19 3/4"	16 1/2"	41"	36 1/2"	25 1/2"	12"	39 lb
53010250	Detachable Full Arm, Elevating Legrest	16"	16"	19 3/4"	16 1/2"	41"	36 1/2"	25 1/2"	12"	39 lb
53010320	Detachable Desk Arm, Swingaway Footrest	20"	16"	19 3/4"	16 1/2"	41"	36 1/2"	29 1/2"	12"	39 lb
53010330	Detachable Desk Arm, Elevating Legrest	20"	16"	19 3/4"	16 1/2"	41"	36 1/2"	29 1/2"	12"	39 lb
53010340	Detachable Full Arm, Swingaway Footrest	20"	16"	19 3/4"	16 1/2"	41"	36 1/2"	29 1/2"	12"	39 lb
53010350	Detachable Full Arm, Elevating Legrest	20"	16"	19 3/4"	16 1/2"	41"	36 1/2"	29 1/2"	12"	39 lb
53010360	Adjustable Height Desk Arm, Swingaway Footrest	20"	16"	19 3/4"	16 1/2"	41"	36 1/2"	29 1/2"	12"	39 lb
53010370	Adjustable Height Desk Arm, Elevating Legrest	20"	16"	19 3/4"	16 1/2"	41"	36 1/2"	29 1/2"	12"	39 lb

NOTE: Everest & Jennings maintains the policy of continual product improvement, therefore, we reserve the right to make changes without notice. All dimensions are ±1/4".
*Without Front Rigging. **With axles and forks in lower position.



Graham-Field Health Products
2935 Northeast Parkway
Atlanta, Georgia 30360
2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com