

HOME DIALYSIS AND INFUSION EQUIPMENT

Patients electing to dialyze from the comfort of their home maintain a better quality of life. Graham-Field® makes a wide variety of equipment that can help make a patient's home dialysis and treatment sessions easier, safer and more comfortable.

LUMEX Recliners

by graham-field

All recliners shown are ideal for home dialysis and infusion treatment.

- Designed and built for long term seating comfort
- Tente® 5 inch (12.70 cm) casters: 3 locking, 1 directional provide ease of maneuverability
- Broad palette of colors and fabrics
- Made in the USA

Deluxe Preferred Care Recliner-Wide

FR565WG Effortlessly clean or replace easily removed quick-release upholstery.

- Infinite recline positioning within a range
- Available in smaller Preferred Care Recliner (FR565G)
- Extra wide armrests
- 24 inch (60.96 cm) wide seat
- 400 pound (181.44 Kg) weight capacity - EVENLY DISTRIBUTED



Clinical Care Recliner

FR577RG Versatile recliner featuring adjustable top headrest section and three positional comfort.

- Two fold-down smooth surface side tables
- Adjustable headrest and wider, softer armrest
- Trendelenburg position
- Available with Heat & Massage (FR577GH)
- 21.50 inch (54.61 cm) wide seat
- 290 pound (131.54 Kg) weight capacity - EVENLY DISTRIBUTED



Shown with optional side table, IV mount and pole

Ortho-Biotic II Clinical Recliner

FR597 Ultimate comfort with infinite recline positioning and independent leg rest operation.

- Self-Storing footplate
- Fully Upholstered (P) or Wood Capped (G) arms and optional Heat and Massage (H) available
- Trendelenburg position
- 23 inch (58.42 cm) wide seat
- 350 pound (158.76 Kg) weight capacity - EVENLY DISTRIBUTED



See the Lumex Seating Swatch Card for the full selection of colors and fabrics.

graham field

HOME DIALYSIS AND INFUSION EQUIPMENT

Patients who perform dialysis at home generally have a better lifestyle. It is important to have a clean and well organized area in which to perform your home dialysis. Graham-Field's patient room equipment can help to organize your "treatment space".



grafco Heating Pad 600 Soothing Dry Heat

4507-1 Many people undergoing their treatments feel cold. The Grafco® Heating Pad is the perfect solution.

- Moist / Dry heat brings comforting relief to sore, aching muscles and joints
- Three automatically controlled heat settings
- Removable, washable outer cover
- Pad Size: 12x15 inches (30.48x 38.10 cm)



LUMISCOPE
by graham-field

Automatic Blood Pressure Monitor

1133 During dialysis, it is important to measure your blood pressure and pulse rate. Lumiscope® offers quality, easy-to-use automatic blood pressure monitors.

- Memory feature - 2 user / 60 memory sets each
- Includes: AC adapter, case and 4 AA batteries
- Irregular heartbeat detection with Pulse Monitor



LUMEX
by graham-field

Overbed Table Non-Tilt Top

GF8902 Lumex OBT is Ideal for patient convenience and placement of treatment prep supplies.

- Walnut, wood-grain laminated 0.75 inch (1.905 cm) particle-board core with plastic edge molding
- Twin-wheel swivel casters: two locking, two non-locking
- Height adjusts infinitely between 28- 44 inches (71.12 - 111.76 cm)



LUMEX
by graham-field

Composite Overbed Tables

- Height adjusts infinitely between 28.3 - 41.7 inches (71.88 - 105.92 cm)
- Moisture impermeable, molded composite top
- Twin-wheel swivel casters: two locking, two non-locking
- Supports up to 50 pounds (22.68 Kg)
- **EVENLY DISTRIBUTED**
- Cleans and disinfects easily


GF890PS Silver Vein Coated Base

GF890P Chrome Plated Base

grafco IV Pole - 4 Leg, 2 Hook

GF7012-1 Grafco IV Pole is perfect for hanging solution bags.

- Swivel ball bearing casters facilitate easy maneuverability on all surfaces
- Height adjusts between 47- 84 inches (119.38 - 213.36 cm)



To see specifications and the complete line of products visit www.grahamfield.com.

1.770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved. Lumex, Grafco and Lumiscope are trademarks of GF Health Products, Inc. GF Health Products, Inc. is an ISO 13485 Certified Company.

graham field

GF16000-42Rev/B19